



# NEW PANTHER SMART STAPLER

Smart reload identification  
Smart tissue recognition  
Smart firing control





One touch firing  
and opening  
system

### SMART HANDLE

Universal handle ready  
to use with the most  
extended portfolio of reloads

### ELECTRONICALLY ASSISTED

Chip incorporated in the reloads  
for the most precise firing



Innovative 45 degrees  
Non-Stop articulation system

Universal handle with  
better design and ergonomics  
for easy 360 degrees knobbing





### Technology P2G

- Double Gripping Technology
- Lateral diffusion of the fluids
- Better tissue compression

Curved tip Reloads  
specially indicated for thin  
tissue and vessels



The P2G technology  
range guarantees  
better stapling delivering  
hemostasis and safety



# NEW GENERATION ENDO STAPLING SYSTEM

Product Name	Product	Description		Color
Smart Stapler	PEAS	Powered Endoscopic Cutter Stapler with Non-Stop articulation system (Short Handle)		
	PEAM	Powered Endoscopic Cutter Stapler with Non-Stop articulation system (Standard Handle)		
	PEAL	Powered Endoscopic Cutter Stapler with Non-Stop articulation system (Extra long Handle)		
Endo Articulating Reloading	CADF30T	PANTHER Endo Articulating reloading unit 30 mm with P2G Technology and 2.5mm staple height (White)	Articulated reload	○ White
	CADF30N	PANTHER Endo Articulating reloading unit 30 mm with P2G Technology and 3.5mm staple height (Blue)	Articulated reload	● Blue
	CADF45T	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 2.5mm staple height (White)	Articulated reload	○ White
	CADF45N	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 3.5mm staple height (Blue)	Articulated reload	● Blue
	CADF-45S	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 4.2mm staple height (Gold)	Articulated reload	● Gold
	CADF-45D	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 4.8mm staple height (Green)	Articulated reload	● Green
	CADF-45R	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 5.0mm staple height (Black)	Articulated reload	● Black
	CADF-60T	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 2.5mm staple height (White)	Articulated reload	○ White
	CADF60N	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 3.5mm staple height (Blue)	Articulated reload	● Blue
	CADF60S	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 4.2mm staple height (Gold)	Articulated reload	● Gold
CADF60D	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 4.8mm staple height (Green)	Articulated reload	● Green	
CADF-60R	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 5.0mm staple height (Black)	Articulated reload	● Black	