

Endo Linear Cutter Stapler


WITH NON-STOP ARTICULATING LEVER

*Innovative 45 degrees
Non-Stop articulation system*

*Ergonomic
handle compatible
with all Panther
reloading units!*

*Loading button on the body
of the instrument*

*Incorporated grasping
mode – grasping the
tissue works by either
pushing or pulling
back the loop handle*



ENDO LINEAR CUTTER STAPLER

ERGONOMIC HANDLE COMPATIBLE WITH ALL PANTHER RELOADING UNITS!

<i>Product</i>	<i>Description</i>
CEAC-30	PANTHER ENDO Linear Cutter Stapler, standard handle with non-stop Articulation system and grasping mode
CEAC-45	PANTHER ENDO Linear Cutter Stapler, short handle with non-stop Articulation system and grasping mode
CEAC-60	PANTHER ENDO Linear Cutter Stapler, extra long handle with non-stop Articulation system and grasping mode

The information presented in this material is intended for product description and promotion only. Panther Healthcare has the relevant power of interpretation. Always refer to the product label and/or instruction for use before using any Panther Healthcare product. For further details and updates please refer to Panther Healthcare website: www.pantherhealthcare.com