

# Endo Linear Cutter Stapler


WITH NON-STOP ARTICULATING LEVER

*Innovative 45 degrees  
Non-Stop articulation system*

*Ergonomic  
handle compatible  
with all Panther  
reloading units!*

*Loading button on the body  
of the instrument*

*Incorporated grasping  
mode – grasping the  
tissue works by either  
pushing or pulling  
back the loop handle*



# ENDO LINEAR CUTTER STAPLER

*ERGONOMIC HANDLE COMPATIBLE WITH ALL PANTHER RELOADING UNITS!*

<i>Product</i>	<i>Description</i>
CEAC-30	PANTHER ENDO Linear Cutter Stapler, standard handle with non-stop Articulation system and grasping mode
CEAC-45	PANTHER ENDO Linear Cutter Stapler, short handle with non-stop Articulation system and grasping mode
CEAC-60	PANTHER ENDO Linear Cutter Stapler, extra long handle with non-stop Articulation system and grasping mode

*The information presented in this material is intended for product description and promotion only. Panther Healthcare has the relevant power of interpretation. Always refer to the product label and/or instruction for use before using any Panther Healthcare product. For further details and updates please refer to Panther Healthcare website: [www.pantherhealthcare.com](http://www.pantherhealthcare.com)*