



Model	Description	Staple Line Length	Cutting Line Length	Staple Height	Staple Formation
CAAB-40	Disposable Linear Cutter Stapler	43mm	40mm	4.8mm	2.0mm
CAZB-40D	Reloading unit for Linear Cutter Stapler	43mm	40mm	4.8mm	2.0mm



Because accessibility matters!!!

Linear Stapler with Knife

Low Anterior Resection

Especially suitable for lower colorectal surgery: it combines horizontal stapling and vertical cutting and supports multiple firing, making operations more flexible.

Because
Accessibility
Matters!




Linear Stapler with Knife



Straight cut edge and staple line
-Well-formed staple line ensuring
smooth blood circulation

Longer shaft makes the device easier
to reach bottom of pelvis
to conduct low rectal surgery

The **smallest head** and
the **largest jaw** aperture
on the market.
Facilitates access and
mobility to place the sta-
pler and the tissue just
exactly where it's
needed



Special design: with its own cutting blade, cutting
and stapling can be done at the same time. The cutting
direction is perpendicular to stapling one; first stapling
then cutting ensures a smoother cut on thicker tissue.

With ergonomic design, it requires **less force on the handle.**