



Panther[®]
Healthcare

A better life through innovative
surgical technology

Panther Tack Hernia Fixation

Panther Tack
a 5mm fixation
device for a secure
laparoscopic
hernia repair



Panther Tack is a safe Hernia
fixation device without
needles or prongs.
It reduces the trauma risks
in Lap operations



Ergonomic design of the handle for
an easy and smooth firing



Panther Tack is presented
in a 30 titanium tack
configuration with a 5mm shaft
and 355mm length

www.pantherhealthcare.com

PRODUCT INFORMATION

Product Code	Description
HRDN0530	Panther Fixation Device with 30 Titanium tacks and 5mm shaft



The information presented in this material is intended for product description and promotion only. Panther Healthcare has the relevant power of interpretation. Always refer to the product label and/or instruction for use before using any Panther Healthcare product. For further details and updates please refer to Panther Healthcare website: www.pantherhealthcare.com